




MORGAN STEER
ORTHOPAEDICS LIMITED

COMMITTED...

TIME HONoured DESIGNS

PROVEN PERFORMANCE



to creating the perfect bowl mix

The Morgan Steer range provides the perfect solution when choosing the right vacuum mixing system.

☑ Quality

☑ Performance

☑ Value

FEATURE	BENEFIT
<ul style="list-style-type: none">Planetary Gearing (Rotational axis paddle).	<ul style="list-style-type: none">Provides complete coverage of the bowl ensuring homogenous, reproducible mix ¹. (Also see document Ref: FB-B-001-APPX).
<ul style="list-style-type: none">Attention to the Highest Quality – Extremely cost effective.	<ul style="list-style-type: none">Offers healthcare providers with quality, effective medical devices at a price point that reflects the challenges and opportunities of the current healthcare environment.
<ul style="list-style-type: none">Closed system during mixing.	<ul style="list-style-type: none">Minimal exposure for staff to monomer fumes in the operating theatre environment.
<ul style="list-style-type: none">Optimum vacuum of 550mmHg.	<ul style="list-style-type: none">Has been shown to give improved mechanical properties, maximising the fatigue strength of the cement ².
<ul style="list-style-type: none">Enhanced Spatula.	<ul style="list-style-type: none">Spatula precisely matches the bowl inner surface profile to help minimise cement wastage.
<ul style="list-style-type: none">Disposable Curette.	<ul style="list-style-type: none">Specifically designed for scraping, shaping and removing excess bone cement.
<ul style="list-style-type: none">Latex Free.	<ul style="list-style-type: none">Protects against latex allergy.
<ul style="list-style-type: none">High volume bowl.	<ul style="list-style-type: none">Allows 40g to 120g of all types of bone cement.
<ul style="list-style-type: none">Charcoal/Microbiological filter.	<ul style="list-style-type: none">Reduces fumes to levels significantly below the HSE guidelines.
<ul style="list-style-type: none">Pictorial Instructions for use.	<ul style="list-style-type: none">Enhanced, very clear step by step guide/ additional training aid for all operating theatre personnel.
<ul style="list-style-type: none">Modern ergonomic design.	<ul style="list-style-type: none">Meeting the current needs of the user but retaining time honoured functional design with a robust emphasis on clinical outcome.
<ul style="list-style-type: none">Individually boxed.	<ul style="list-style-type: none">Ensures no compromising of the inner and outer peel pouches. Meets the rigorous, current CE requirements of packaging and transportation for products in this medical device class.
<ul style="list-style-type: none">Next day delivery.	<ul style="list-style-type: none">Meeting purchasing and storage needs of healthcare providers.
<ul style="list-style-type: none">NHS Supply Chain.	<ul style="list-style-type: none">Alternative purchasing route for healthcare professionals offering choice, value and flexibility.